THERMOMETER - THERMOMÈTRE

INVENTUM

SINCE 1908

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ENGLISH

Legislation requires that all electrical and electronic equipment must be collected for reuse and recycling. Electrical and electronic equipment marked with the symbol indicating separate collection of such equipment must be returned to a municipal waste collection point.


Device classification type BF


Please read the instructions for use!

DESCRIPTION OF THE DEVICE	PAGE 25
1. SAFETY INSTRUCTIONS	PAGE 26
2. INTRODUCTION	PAGE 27
3. TAKING THE TEMPERATURE	PAGE 28
4. MEMORY	PAGE 28
5. CHANGING THE MEASUREMENT UNIT	PAGE 29
6. EXPLANATION OF SYMBOLS/ERROR MESSAGES	PAGE 29
7. CHANGING THE BATTERY	PAGE 30
8. CLEANING, STORAGE & DISPOSAL	PAGE 30
9. TECHNICAL DATA	PAGE 31
GARANTIE & SERVICE	PAGE 35

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

front side

rear side

1. Sensor tip
2. LCD-Display
3. On/ Off/ Measurement button
4. Protective cap
5. Memory button **M**
6. Battery compartment cover
7. Battery compartment unlocking

1 | SAFETY INSTRUCTIONS

- **Read the instruction manual carefully and completely before using the appliance and keep them for later use.**
- This device is a sensitive electronic device. Please treat it carefully and do not expose it to any mechanical impacts.
- Do not expose the thermometer to direct sunlight.
- The device should be in the room in which the measurement will be made for at least 30 minutes before use.
- The thermometer, with the exception of the measuring tip, is NOT waterproof. For this reason, direct contact with water or other liquids should be avoided.
- Clean the waterproof measuring tip after each use with a soft cloth moistened with disinfectant.
- Check before each use that the lens in the measuring tip is intact. If it is damaged, please contact the dealer or service address.
- The thermometer was constructed for practical use but cannot replace a visit to the doctor.
- This device is not intended for commercial or clinical use.
- There are people in whom different measurement results are obtained in the left and right ear. In order to record temperature changes, always measure in the same ear of the same person.
- If you have been lying on one ear for some time, the temperature is slightly raised. Wait a little while or measure in the other ear.
- As ear wax can affect the measurement, you should clean the ear before measuring if necessary.
- The ear thermometer may be used by children only under adult supervision. Measurement is usually possible from the age of 6 months and up.
In infants under 6 months, the ear passage is still very narrow so that the temperature of the eardrum can often not be recorded and the measurement result displayed is too low.
- The measurement must not be taken in an ear affected by inflammatory diseases (e.g. discharging pus or secretion), after possible ear injuries (e.g. eardrum damage) or in the healing period after operative procedures. In all of these cases, please talk to your doctor.
- The device is intended only for the purpose given in these directions for use.
- Repairs must be carried out only by authorised service agents. Otherwise the guarantee becomes void.
- **Never immerse the thermometer in water or any other liquid.**
- This device complies with EU Directive 93/42 EEC.
- This device complies with standard EN 12470-5 Medical Thermometers, Requirements for Infrared Ear Thermometers (with maximum device).
- If you have further questions about using our devices, please contact your dealer.

2 | INTRODUCTION

What you should know about using the Inventum ear thermometer.

The thermometer is suitable for measuring body temperature in the ear. The temperature varies depending on the part of the body where the measurement is taken. The difference in healthy persons between different parts of the body can be between 0.2-1 °C.

The normal temperature ranges measured with an ear thermometer are between 36,0 to 37,8 °C. That means that the temperature measured in the ear reflects the actual body temperature. This has the advantage that temperature changes are identified faster when measured in the ear than with rectal measurement. This is particularly important when the patient's temperature is rising or falling.

Tip: Temperatures measured with different thermometers cannot be compared with one another. You should therefore tell your doctor or bear in mind if you are diagnosing yourself what thermometer you used to take your temperature and in what part of the body.

The temperature of a healthy person is also influenced by the following factors:

- The person's individual metabolism
- Age (body temperature is higher in babies and small children and falls with increasing age. Greater temperature fluctuations occur faster and more often with children.)
- Clothing
- The outside temperature
- The time of day (body temperature is lower in the morning and increases during the day towards evening)
- Preceding physical and, to a lesser extent, mental activity.

Tip: Taking the temperature gives a measurement that provides information about a person's current body temperature. If you are uncertain about interpreting the results or if the result is abnormal (e.g. fever), you should consult your doctor.

This also applies in the case of slight temperature changes if there are other symptoms of illness such as agitation, severe sweating, flushed skin, fast pulse rate, tendency to collapse etc.

3 | TAKING THE TEMPERATURE

This thermometer is designed for use in the ear (ear passage / eardrum) usually for persons older than 6 months. With children under 6 months, the ear passage is still very narrow so that the temperature of the eardrum can often not be recorded and the measurement result displayed is too low. Before taking the temperature, remove the protective cap and ensure that the sensor tip and also the ear passage are clean. If there are signs of acute inflammation (discharge of pus or secretion, pain), in the case of injuries or directly after operations on the ear, the non-affected ear should be chosen. Otherwise, the measurement result can be incorrect. Use of the thermometer by different persons can be inappropriate in certain acute infectious diseases because of the possible spread of germs despite cleaning and disinfection. Discuss this with your doctor in the individual case. The thermometer must only be used without disposable covers. There is already a battery in it already so it can be used immediately when operated for the first time.

- Switch on the thermometer with the 'ⓘ' button. After a short self-test and a short beep the thermometer is ready for measuring the temperature in the ear.
- Because the ear canal is slightly curved, pull the ear back and up slightly before inserting the sensor tip.

This is important to allow the sensor tip to be pointed directly at the ear drum.

- Insert the sensor tip carefully and press the button 'ⓘ' for about 1 sec.
- The end of the measuring time is signalled by a short beep.
- Now you can read the measured temperature.
- The smiley is also illuminated according to the test result:

😊 indicates that the body temperature is in the normal range,
☹️ is displayed at a temperature of or higher than 38°C, i.e. fever alarm.

- The time between measurements should be at least 10 seconds. This ensures that the measurement is more accurate.

4 | MEMORY

The temperature last measured, is saved automatically. There are 10 memory buttons available for this. The last temperature measurements can be recalled by pressing the 'M' button on the bottom of the thermometer.

5 | CHANGING THE MEASUREMENT UNIT

The temperature can be displayed in degrees Celsius ($^{\circ}\text{C}$) or degrees Fahrenheit ($^{\circ}\text{F}$). Press the button ‘ C/F ’ to change the measuring unit. Without releasing it press the ‘ M ’ button to switch to the desired measuring unit. All saved values are displayed in the new measuring unit.


6 | EXPLANATION OF SYMBOLS/ ERROR MESSAGES

Message d'erreur	Problème	Solution
<i>Er</i>	Ambient temperature below 16°C or above 35°C ($< 60,8^{\circ}\text{F}$, $> 95^{\circ}\text{F}$)	The ambient temperature must be between 16°C and 35°C ($60,8^{\circ}\text{F}$, 95°F).
<i>Hi</i>	The temperature measured is higher than 43°C ($109,4^{\circ}\text{F}$).	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your dealer or customer service.
<i>Lo</i>	The temperature measured is lower than 34°C ($93,2^{\circ}\text{F}$).	Operate the thermometer only between the specified temperature ranges. If necessary, clean the sensor tip. In the event of a repeated error message, contact your dealer or customer service.

7 | CHANGING THE BATTERY

The type 3V CR2032 is already fitted and depending on usage has a life of around 3000 measurements. When the battery gets weak, the battery warning symbol appears . It is still possible to measure the temperature. Battery must be replaced. If the thermometer displays the battery symbol  and 'Lo', the battery must be replaced.

- Remove the screw from the battery cover.
- Open the battery cover and remove the original battery.
- Take a new CR2032 lithium battery and align it with the positive terminal (anode) facing upwards.
- Place the battery into the battery compartment.
- Replace the battery cover back tightly.
- Screw the battery cover back tightly.

Used batteries should not go into domestic refuse. You are legally obliged to dispose of the batteries correctly. Dispose of them through your electrical dealer or local recycling centre. **Note:** you will find these symbols on batteries containing toxic substances: Pb = battery contains lead, Cd = battery contains cadmium, Hg = battery contains mercury. The battery in this device does not contain toxic substances.

8 | CLEANING, STORAGE & DISPOSAL

Clean the waterproof sensor tip after every use. Use a clean cloth or cotton bud which is moistened with disinfectant, alcohol or warm water. To clean the entire device, use a soft cloth slightly moistened with a mild soapy solution.

Water must never get into the device. If water gets into the device, remove the battery immediately. Only use the device again when it's completely dry.

Do not use any harsh cleaning products. Never immerse the device in water. Only the tip is waterproof.

Always store the thermometer with the cap on in order to protect the sensor. When putting the cap on the device, align the arrow on the cap with the arrow on the housing.

The device must not be stored or used at too high or low a temperature or humidity (see technical specifications), in sunlight, in association with an electrical current or in dusty locations. Otherwise, inaccuracies can occur.

If prolonged storage is intended, you should remove the battery.

10 | TECHNICAL DATA

Model	MT22
Measurement range	34 °C - 43 °C (93,2 °F - 109,4 °F)
Laboratory measurement precision	±0,2 °C (±0,4 °F) between 35.5 °C - 42 °C (95.9 °F - 107.6 °F) ±0,3 °C (±0,5 °F) below 35.5 °C and above 42 °C
Clinical repeat precision	0,31 (children), 0,28 (adults)
Measurement duration	1 to 2 seconds
Interval between 2 measurements	At least 10 seconds
Measurement units	°Celsius (°C) or °Fahrenheit (°F)
Operating conditions	16 °C to 35 °C (60.8 °F - 95 °F) with a relative humidity of up to 85% (non-condensing)
Storage conditions	-25 °C to 55 °C (-13 °F - 131 °F) with a relative humidity of up to 85% (non-condensing)
Dimensions	54 x 44 x 158 mm
Weight	65.5 g including battery
Battery	1 lithium battery (type 3V CR-2032)
Memory	For 10 measurements
Explanation of symbols	Device classification type BF 
	Please read the instructions for use! 

Note: When the device is used outside the specifications, perfect functioning cannot be guaranteed!

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GUARANTEE & AFTER SALES SERVICE

- This product is guaranteed for a period of 2 years from the date of purchase against mechanical and electrical defect.
- Guarantee can only be given if proof of purchase (e.g. till receipt) can be shown. Therefore keep it in a safe place. The guarantee period starts at date of purchase.
- Guarantee is not valid in case of:
 - normal wear;
 - damage through misuse use or when the instructions were not followed correctly;
 - lack of maintenance;
 - not following the operation- and safety instructions;
 - incompetent montage or repairs made by yourself or others;
 - using accessories and spare parts which are not suitable for the appliance; for assessment by our technical department.
- The guarantee is not valid for:
 - normal wear;
 - batteries;
 - (interchangeable) signal- and controllamps for assessment by our technical department.

Return to your supplier or send the appliance, including a copy proof of purchase and description of the complaint, stamped to:

Inventum Holland BV
Generatorstraat 17
3903 LH VEENENDAAL
THE NETHERLANDS

new E.T.P. bvba
Passtraat 121
9100 Sint Niklaas
BELGIUM

- Pack it carefully (preferably in the original carton). Ensure unit is clean.
- Enclose your name and address and quote model number on all correspondence.

When the guarantee is no longer valid or when the period of guarantee has passed, we will charge any costs that flow out of an assessment and/or repairing.

On request we will supply a statement before repairing. Without requesting a statement, we assume that you, by sending the appliance, automatically order to repair the appliance.

A guarantee period of two months is valid for repairs done by our technical department and the used spare parts, counting from the repairedate. When a malfunction occurs within this period which directly flows out of the repairs and the used spare parts, the consumer has the right of repairs without any costs.

Inventum Holland BV
Generatorstraat 17, 3903 LH , Veenendaal, Nederland
Tel.: 0318-542222, Fax.: 0318-542022, Mail: helpdesk@inventum.eu

Inventum Belgium BVBA
Lindenstraat 32, 3140, Keerbergen, België
Tel: 015-512653, Fax.: 015-529718, Mail: info@inventumbelgium.be

Website: www.inventum.eu